

Chikungunya IgM/IgG Rapid Test Cassette (Serum / Plasma/Whole Blood)

Cat no. CHIK-441

INTENDED USE

The Chikungunya IgM/IgG Rapid Test cassette is a lateral flow chromatographic immunoassay for the qualitative detection of IgM/IgG anti-chikungunya virus (CHIK) in human serum, plasma or whole blood. It is intended to be used as a screening test and as an aid in the diagnosis of infection with CHIK. Any reactive specimen with the Chikungunya IgM/IgG Rapid Test must be confirmed with alternative testing method(s) and clinical findings.

SUMMARY AND EXPLANATION OF THE TEST

Chikungunya is a rare viral infection transmitted by the bite of an infected *Aedes aegypti* mosquito. It is characterized by a rash, fever, and severe joint pain (arthralgias) that usually lasts for three to seven days. The name is derived from the Makonde word meaning "that which bends up" in reference to the stooped posture developed as a result of the arthritic symptoms of the disease. It occurs during the rainy season in tropical areas of the world, primarily in Africa, South-East Asia, southern India and Pakistan^{1,2}.


The symptoms are most often clinically indistinguishable from those observed in dengue fever. Indeed, dual infection of dengue and chikungunya has been reported in India³. Unlike dengue, hemorrhagic manifestations are relatively rare and most often the disease is a self limiting febrile illness. Therefore it is very important to clinically distinguish dengue from CHIK infection.

CHIK is diagnosed based on serological analysis and viral isolation in mice or tissue culture. An IgM immunoassay is the most practical lab test method⁴.

The Chikungunya IgM/IgG Rapid Test utilizes recombinant antigens derived from its structure protein⁵, it detects IgM/IgG anti-CHIK in patient serum or plasma within 15 minutes. The test can be performed by untrained or minimally skilled personnel, without cumbersome laboratory equipment.

TEST PRINCIPLE

The Chikungunya IgM/IgG Rapid Test is a lateral flow chromatographic immunoassay. The test cassette consists of: 1) a burgundy colored conjugate pad containing recombinant antigen conjugated with colloidal gold (Chikungunya conjugates) and rabbit IgG-gold conjugates, 2) a nitrocellulose membrane strip containing two test bands (M and G bands) and a control band (C band). The M band is pre-coated with monoclonal anti-human IgM for the detection of IgM anti-Chikungunya, G band is pre-coated with reagents for the detection of IgG anti-Chikungunya, and the C band is pre-coated with goat anti rabbit IgG.



When an adequate volume of test specimen is dispensed into the sample well of the test cassette, the specimen migrates by capillary action across the cassette. Anti-CHIK IgM if present in the specimen will bind to the CHIK conjugates. The immunocomplex is then captured on the membrane by the pre-coated anti-human IgM antibody, forming a burgundy colored M band, indicating a CHIK IgM positive test result.

Anti-CHIK IgG if present in the specimen will bind to the CHIK conjugates. The immunocomplex is then captured by the pre-coated reagents on the membrane, forming a burgundy colored G band, indicating a CHIK IgG positive test result.

Absence of any test bands (M and G) suggests a negative result. The test contains an internal control (C band) which should exhibit a burgundy colored band of the immunocomplex of goat anti rabbit IgG/rabbit IgG-gold conjugate regardless of the color development on any of the test bands. Otherwise, the test result is invalid and the specimen must be retested with another device.

REAGENTS AND MATERIALS PROVIDED

- Each sealed in a foil pouch with three items inside:
 - One cassette device.
 - One pipette dropper.
 - One desiccant.
- Sample Diluent (1 vial, 4 mL)
- One package insert (instruction for use).

MATERIALS REQUIRED BUT NOT PROVIDED

- Clock or Timer

WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use

- This package insert must be read completely before performing the test. Failure to follow the insert gives inaccurate test results.
- Do not open the sealed pouch, unless ready to conduct the assay.
- Do not use expired devices.
- Bring all reagents to room temperature (15°C-30°C) before use.
- Do not use the components in any other type of test kit as a substitute for the components in this kit.
- Do not use hemolyzed blood specimen for testing.
- Wear protective clothing and disposable gloves while handling the kit reagents and clinical specimens. Wash hands thoroughly after performing the test.
- Users of this test should follow the US CDC Universal Precautions for prevention of transmission of HIV, HBV and other blood-borne pathogens.

- Do not smoke, drink, or eat in areas where specimens or kit reagents are being handled.
- Dispose of all specimens and materials used to perform the test as biohazardous waste.
- The testing results should be read within 15 minutes after a specimen is applied to the sample well or sample pad of the device. Read result after 15 minutes may give erroneous results.
- Do not perform the test in a room with strong air flow, ie. an electric fan or strong air-conditioning.

REAGENT PREPARATION AND STORAGE INSTRUCTIONS

All reagents are ready to use as supplied. Store unused test device unopened at 2°C-30°C. If stored at 2°C-8°C, ensure that the test device is brought to room temperature before opening. The test device is stable through the expiration date printed on the sealed pouch. Do not freeze the kit or expose the kit over 30°C.

SPECIMEN COLLECTION AND HANDLING

Consider any materials of human origin as infectious and handle them using standard biosafety procedures.

Plasma

- Collect blood specimen into a lavender, blue or green top collection tube (containing EDTA, citrate or heparin, respectively in Vacutainer®) by veinpuncture.
- Separate the plasma by centrifugation.
- Carefully withdraw the plasma into new pre-labeled tube.

Serum

- Collect blood specimen into a red top collection tube (containing no anticoagulants in Vacutainer®) by veinpuncture.
- Allow the blood to clot.
- Separate the serum by centrifugation.
- Carefully withdraw the serum into a new pre-labeled tube.

Test specimens as soon as possible after collecting. Store specimens at 2°C-8°C if not tested immediately.

Store specimens at 2°C-8°C up to 5 days. The specimens should be frozen at -20°C for longer storage.

Avoid multiple freeze-thaw cycles. Prior to testing, bring frozen specimens to room temperature slowly and mix gently. Specimens containing visible particulate matter should be clarified by centrifugation before testing.

Do not use samples demonstrating gross lipemia, gross hemolysis or turbidity in order to avoid interference on result interpretation.

Blood

Drops of whole blood can be obtained by either finger tip puncture or veinpuncture. Do not use any hemolyzed blood for testing.

Whole blood specimens should be stored in refrigeration (2°C-8°C) if not tested immediately. The specimens must be tested within 24 hours of collection.

ASSAY PROCEDURE

Step 1: Bring the specimen and test components to room temperature if refrigerated or frozen. Mix the specimen well prior to assay once thawed.


Step 2: When ready to test, open the pouch at the notch and remove device. Place the test device on a clean, flat surface.

Step 3: Be sure to label the device with specimen's ID number.

For whole blood test

Apply 1 drop of whole blood (about 40-50 µL) into the sample well.

Then add 1 drop (about 35-50 µL) of Sample Diluent immediately.



For serum or plasma test

Fill the pipette dropper with the specimen.

Holding the dropper vertically, dispense 1 drop (about 30-45 µL) of specimen into the sample well making sure that there are no air bubbles.

Then add 1 drop (about 35-50 µL) of Sample Diluent immediately.



Chikungunya IgM/IgG Rapid Test Cassette (Serum / Plasma/Whole Blood)

Cat no. CHIK 441

Step 5: Set up timer.

Step 6: Results can be read in 15 minutes. Positive results can be visible in as short as 1 minute.


Don't read result after 15 minutes. To avoid confusion, discard the test device after interpreting the result.

QUALITY CONTROL

Although the testing device contains an internal quality control (pink colored band in the control region), good laboratory practice recommends the daily use of an outside control to ensure proper testing device performance. Quality control samples should be tested according to the standard quality control requirements established by your laboratory.

INTERPRETATION OF ASSAY RESULT

- 1. NEGATIVE RESULT:** If only the C band is present, the absence of any burgundy color in the both test bands (M and G) indicates that no anti-CHIK antibody is detected. The result is negative.



- 2. POSITIVE RESULT:**

- 2.1** In addition to the presence of C band, if only M band is developed, the test indicates for the presence of anti-CHIK IgM. The result is positive.



- 2.2** In addition to the presence of C band, if only G band is developed, the test indicates for the presence of anti-CHIK IgG. The result is positive.




- 2.3** In addition to the presence of C band, both M and G bands are developed, the test indicates for the presence of anti-CHIK IgG and IgM. The result is also positive.



Samples with positive results should be confirmed with alternative testing method(s) and clinical findings before a positive determination is made.

- 3. INVALID:** If no C band is developed, the assay is invalid regardless of any burgundy color in the test bands as indicated below. Repeat the assay with a new device.



LIMITATIONS OF TEST

- The Assay Procedure and the Assay Result Interpretation must be followed closely when testing the presence of IgM anti-CHIK in serum, plasma or whole blood from individual subjects. Failure to follow the procedure may give inaccurate results.
- The Chikungunya IgM/IgG Rapid Test is limited to the qualitative detection of IgM/IgG anti-CHIK in human serum, plasma or whole blood. The intensity of the test band does not have linear correlation with the antibody titer in the specimen.
- A negative result for an individual subject indicates absence of detectable IgM/IgG anti-CHIK. However, a negative test result does not preclude the possibility of exposure to or infection with CHIK.
- A negative result can occur if the quantity of IgM/IgG anti-CHIK present in the specimen is below the detection limits of the assay, or the antibodies that are detected are not present during the stage of disease in which a sample is collected.
- Some specimens containing unusually high titer of heterophile antibodies or rheumatoid factor may affect expected results.
- The results obtained with this test should only be interpreted in conjunction with other diagnostic procedures and clinical findings.

REFERENCES

- Shah KV, Gibbs CJ Jr, Banerjee G. Virological investigation of the epidemic of haemorrhagic fever in Calcutta: isolation of three strains of Chikungunya virus. Indian J Med Res 1964; 52 :676-83.
- Powers AM, Brault AC, Tesh RB, Weaver SC. Re-emergence of Chikungunya and O'nyong-nyong viruses: evidence for distinct geographical lineages and distant evolutionary relationships. J Gen Virol 2000;81:471-9
- Myers RM and Carey DE. Concurrent isolation from patient of two arboviruses, Chikungunya and dengue type 2. Science 1967;157:1307-8.
- Thein S, La Linn M, Aaskov J, Aung MM, Aye M, Zaw A, Myint A. Development of a simple indirect enzyme-linked immunosorbent assay for the detection of immunoglobulin M antibody in serum from patients Following an outbreak of chikungunya virus infection in Yangon, Myanmar. Trans R Soc Trop Med Hyg. 1992 86:438-42.
- Yamamoto K, Hashimoto K, Ogata T Structural proteins of Chikungunya virus. Simizu B, J Virol. 1984 Jul;51(1):254-8

MANUFACTURER

ATLAS LINK TECHNOLOGY CO., LTD

Road Xing Min, Guan South Industry Zone,
065500 Langfang City, Hebei Province, CHINA
WEB: <https://www.invitro-test.com>

REVISION:20170316BJ/01/Chik-1