




One Step Strip Style FSH Urine Test

INTRODUCTION

One Step Strip Style FSH Urine Test is a test kit for the determination of FSH (human Follicular Stimulating Hormone) in urine specimens. This test kit is used to obtain a visual result.


TEST PROCEDURE

1. When you are ready to begin testing, open the sealed pouch by tearing along the notch. Remove the test from the pouch.
2. Immerse the strip into the urine with the arrow end pointing towards the urine. Do not immerse past the MAX (maximum) line. Take the strip out after 3 seconds and lay the strip flat on a clean, dry, non-absorbent surface (e.g., mouth of the urine container).
3. Wait for colored bands to appear. Depending on the concentration of FSH in the test specimen, strong positive results may be observed in as short as 40 seconds. However, to confirm negative results, the complete reaction time (5 minutes) is required. Do not read results after 30 minutes.



INTERPRETATION OF RESULTS

- **Negative:** Only one color band appears on the control region or the test band is lighter than the control band. The FSH concentration of the sample is below the detection sensitivity of 25mIU/ml.
- **Positive:** If two color bands are visible and the test band is equal to or darker than the control band. The FSH is present in the sample at or above the detection sensitivity of 25mIU/ml.
- **Invalid:** No visible band at all. Repeat test with a new test kit.



STORAGE AND STABILITY

The test kits can be stored at room temperature (18 to 30°C) in the sealed pouch to the date of expiration. The test kits should be kept away from direct sunlight, moisture, and heat.

PRECAUTION

1. For in vitro diagnostic use only.
2. Do not use test kit beyond expiry date.
3. The test device should not be reused.