


# One Step Midstream Style FSH Urine Test




## INTRODUCTION

One Step Midstream Style FSH Urine Test is a test kit for the determination of FSH (human Follicular Stimulating Hormone) concentration in urine specimens. Test results are read visually without any instrument.



## TEST PROCEDURE


1. When you are ready to begin testing, open the sealed pouch by tearing along the notch. Remove test kit from the pouch.
2. Hold the round end of cover with one hand. Use the other hand to pull out the test device and expose the absorbent.
3. Point the absorbent tip downward; place the absorbent tip in urine stream for at least 10 seconds to be thoroughly wet. Otherwise, you can collect your urine into a clean cup and dip half of the absorbent pad into the urine for at least 10 seconds.



4. Re-cap the device and wait for color bands to appear. Depending on the concentration of FSH, positive result may be observed in as short as 40 seconds. However, to confirm negative results, the complete reaction time (10 minutes) is required. Do not read results after 30 minutes.

## INTERPRETATION OF RESULTS

- **Negative:** Only one color band appears on the control (C) region. Or the test (T) band is lighter than the control band. The FSH concentration of the sample is below 25mIU/ml.
- **Positive:** Distinct color bands appear on the control (C) and test (T) regions. Test line is equal to or darker than the control band. The FSH is present in the sample at or above 25 mIU/ml.
- **Invalid:** No visible band at all. Repeat test with a new test kit.



## STORAGE AND STABILITY

The test kit can be stored at room temperature (4 to 30 °C) in the sealed pouch to the date of expiration. The test kits should be kept away from direct sunlight, moisture and heat.

Manufacturer: Atlas Link Technology Co., Ltd

Web: <https://www.invitro-test.com>