

Atlas Link (Beijing) Technology Co., Ltd

One Step Hepatitis B Surface Antigen (HBsAg) Rapid Test

Instructions For Use

Format: Strip

Specimen: Serum/Plasma

Catalog Number: HBV211

CE 2265

INTENDED USE

One Step Hepatitis B Surface Antigen (HBsAg) Rapid Test is a rapid and convenient immunochromatographic assay for qualitative detection of HBsAg in human serum or plasma sample. It is intended for professional use as an aid in the diagnosis of Hepatitis B virus (HBV) infection. This assay provides only a preliminary result. Clinical expertise and professional judgment should be sought to further evaluate the result of the test

SUMMARY AND PRINCIPLE OF THE ASSAY

Hepatitis B viruses (HBV) have partially double-stranded DNA that causes acute or chronic hepatitis. The viral particles are transmitted through exposure of infectious body fluids or blood, blood transfusion and use of contaminated needles or syringes. Chronic hepatitis may progress to severe outcomes, including cirrhosis and liver cancer (hepatocellular carcinoma). The virus is divided into four major serotypes (adr, adw, ayr, ayw) based on antigenic epitopes presented on its envelope proteins.

Hepatitis B surface antigen (HBsAg) is the first marker to appear in the blood in acute hepatitis B., It can be detected 1 week to 2 months after exposure and 2 weeks to 2 months before the onset of symptoms. Three weeks after the onset of acute hepatitis almost half of the patients will still be positive for HBsAg. In the chronic carrier state, the HBsAg persists for long periods with no seroconversion to the corresponding antibodies. The most commonly used diagnostic and blood screening markers sought is HBsAg. An individual positive for HBsAg is considered to be infected with HBV and is therefore potentially infectious.

One Step HBsAg test is an antigen-capture immunochromatographic assay, it detects the presence of HBsAg in serum or plasma samples. Monoclonal antibodies specifically against HBsAg are 1) conjugated with colloidal gold and deposited on the conjugate pad and 2) immobilized on the test line on the nitrocellulose membrane. When serum or plasma sample is added, it rehydrates the gold-antibody conjugate. If the HBsAg is present in the sample, it will interact with the gold conjugated antibodies. The antigen-antibody-gold complex will migrate towards the test window until the Test Zone (T) where they will be captured by immobilized antibodies, forming a visible red line (Test band) indicating a positive result. If HBsAg are absent in the sample, no red line will appear in the Test Zone (T) indicating a negative result.

To serve as an internal process control, a control line should always appear at Control Zone (C) after the test is completed. Absence of a colored control line in the Control Zone is an indication of an invalid result.

PACKAGE CONTENTS

- Pouch contents: Test Strip, desiccant
- Test instruction

OTHER REQUIRED MATERIALS (NOT PROVIDED)

- Gloves
- Clock or timer

WARNINGS AND PRECAUTIONS

- For professional *in vitro* diagnostic use only.
- Do not reuse.
- Do not use if the pouch seal or its packaging is compromised.
- Do not use after the expiration date shown on the pouch.
- Do not mix and interchange different specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection while handling potentially infectious materials or performing the assay.
- Wash hands thoroughly after finishing the tests.
- Do not eat, drink or smoke in the area where the specimens or kits are being handled.

- Clean up spills thoroughly with appropriate disinfectants.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing procedures.
- Dispose all specimens and used devices in a proper bio-hazard container. The handling and disposal of the hazardous materials should follow local, national or regional regulations.
- Keep out of children's reach.

SPECIMEN PREPARATION

- For serum samples, collect blood in a tube without anticoagulant and allow it to clot.
- For plasma samples, collect blood in a tube containing anticoagulant.
- Separate serum or plasma from blood as soon as possible to avoid hemolysis. Use only clear, non-hemolyzed specimens.
- Testing should be performed immediately after the specimens have been collected. Do not leave the specimens at room temperature for prolonged periods.
- The blood may be stored at 2°C to 8°C for up to three days if the tests cannot be performed immediately. Ensure that the blood samples should be allowed to attain room temperature prior to use.

TEST PROCEDURES

1

Remove the testing device from the foil pouch by tearing at the notch. Hold the strip at the colored end. (Do not touch the arrow end; do not touch test window, the middle part of the strip)



2

Holding the strip vertically, immerse the strip into the specimen with the arrow end pointing towards the specimen. Do not immerse past the MAX line.



3

Take the strip out after a minimum of 10 sec. Lay the strip (MAX side facing up) flat on a clean, dry, nonabsorbent surface.



4


Read the result in 10 minutes. Read results as shown under interpretation of Results.

NOTE: Specimens with high concentrations of HBsAg may produce positive result in as little as 1 minute. Confirm negatives in 20 minutes.



**DO NOT INTERPRET RESULTS
AFTER 30 MINUTES**

RESULT INTERPRETATIONS



Negative

A red colored band appears only at the control region (C), indicating a negative result for HBV infections.

Positive

A clear red control band (C) and a detectable test band (T) appear, indicating a positive result for HBV infections.

Invalid

No visible band appears at the control region (C). Repeat with a new test device. If test still fails, please contact the distributor with the lot number

QUALITY CONTROL

Although the testing device contains an internal quality control (red colored band in the control region), good laboratory practice recommends the daily use of an outside control to ensure proper testing device performance. Quality control samples should be tested according to the standard quality control requirements established by your laboratory.

STORAGE AND STABILITY

- Test device in the sealed pouch should be stored at 2-30°C. Do not freeze the test device.
- The test device should be kept away from direct sunlight, moisture and heat.

LIMITATIONS

- This product is an *in vitro* diagnostic test designed for professional use only.
- Humidity and temperature can adversely affect results.
- The instructions for the use of the test should be followed during testing procedures.
- There is always a possibility that false results will occur due to the presence of interfering substances in the specimen or factors beyond the control of the manufacturer, such as technical or procedural errors associated with the testing.
- Although the test demonstrates superior accuracy in detecting HBsAg, a low incidence of false results can occur. Therefore, other clinically available tests are required in case of questionable results. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

PERFORMANCE CHARACTERISTICS

1. Analytic Sensitivity

One Step HBsAg Test detects HBsAg for major serotypes (adr, adw, ayr, ayw) at concentrations of 1.0 ng/ml.

2. Analytic Specificity

The effect of seromarkers associated with unrelated medical conditions on the specificity of the one Step HBsAg Test was assessed using a panel of specimens. The seromarkers studied were: HCV samples, Dengue samples, the haemolytic samples, rheumatoid factors-contained samples, HIV samples, anti-CMV positive, anti-EBV positive, anti-HSV positive, anti-HAV positive, non-viral liver diseases (non-alcoholic steatohepatitis, fatty liver, drug-induced hepatotoxicity), rubella antibody positive, syphilis serology positive, anti-nuclear antibody positive, C-reactive protein (CRP) positive, antistreptolysin O titre (ASOT) positive. The results demonstrated that One Step HBsAg Test kits have no cross-reactivity with these specimens.

3. Diagnostic Sensitivity and Specificity

Total 554 positive samples were tested 553 positive by HBsAg Test and CE marked Wantai HBsAg EIA; the one (No.216 from ChongYi) tested negative by both HBs Ag Rapid Test and Wantai HBsAg EIA, was further confirmed positive by PCR. The serotype of the false negative sample was unknown. The diagnostic sensitivity of HBsAg Rapid Test was 99.82% (553/554) and could identify the 4 HBV serotypes.

Table 1 Summary of Diagnostic Sensitivity of HBsAg Rapid Test

Serotype	Results of HBs Antigen Rapid Test		Results of CE Marked Wantai HBsAg EIA		Subtotal
	Positive	Negative	Positive	Negative	
adw2	118	0	118	0	118
ayw1	33	0	33	0	33
adr	5	0	5	0	5
ayw2	3	0	3	0	3
ayr	1	0	1	0	1
ayw3	1	0	1	0	1
Unknown serotype	392	1	392	1	393
Subtotal	553	1	553	1	554

Of the 1752 samples, 1750 were tested negative, 2 tested positive (both from blood donors) by HBsAg Rapid Test. All the results which showed inconsistent with CE marked Wantai HBsAg EIA were further confirmed negative by PCR. The diagnostic sensitivity of HBsAg Rapid Test was 99.89% (1750/1752), false positive rate is 0.11%.

Table 2 Summary of Diagnostic Specificity of HBsAg Rapid Test

	Results of HBsAg Rapid Test		Results of CE Marked Wantai HBsAg EIA		Subtotal
	Negative	Positive	Negative	Positive	
Blood Donors	1114	2	1113	3	1116
200 clinical specimens	200	0	200	0	200
200 Pregnant Women	200	0	200	0	200
236 potentially interfering samples	236	0	236	0	236
Subtotal	1750	2	1749	3	1752

4. Interference Study

The following substances and conditions were found not to interfere with the test. List of potentially interfering chemical analytes and concentrations tested are as follow:

Interfering compounds or factors		With Interfering compound added			Without interfering compound		
Name	Spiked final Conc.	Negative serum	HBsAg positive serum (1ng/ml)	HbsAg positive serum (5ng/ml)	Negative serum	HBsAg positive serum (1ng/ml)	HbsAg positive serum (5ng/ml)
Chemical analytes							
Acetaminophen	200 ug/ml	-	+	+	-	+	+
Acetylsalicylic Acid	200 ug/ml	-	+	+	-	+	+
Amikacin	200 ug/ml	-	+	+	-	+	+
Ascorbic acid	200 ug/ml	-	+	+	-	+	+
Aspartame	200 ug/ml	-	+	+	-	+	+
Atropine Sulfate	200 ug/ml	-	+	+	-	+	+
Benzoic Acid	200 ug/ml	-	+	+	-	+	+
Caffeine	200 ug/ml	-	+	+	-	+	+
Deoxyephedrine	200 ug/ml	-	+	+	-	+	+
Dextromethorphan	200 ug/ml	-	+	+	-	+	+
EDTA	800 ug/ml	-	+	+	-	+	+
Gentesic acid	200 ug/ml	-	+	+	-	+	+
Histamine	200 ug/ml	-	+	+	-	+	+
Methaqualone	200 ug/ml	-	+	+	-	+	+
Pendimethazine	200 ug/ml	-	+	+	-	+	+
Penicillin G	200 ug/ml	-	+	+	-	+	+
Quinine	200 ug/ml	-	+	+	-	+	+
Ranitidine	200 ug/ml	-	+	+	-	+	+
Sodium Salicylate	200 ug/ml	-	+	+	-	+	+
Tryptophan	200 ug/ml	-	+	+	-	+	+
Tetracycline	200 ug/ml	-	+	+	-	+	+
Tetrahydrozoline	200 ug/ml	-	+	+	-	+	+
Ethanol	1%	-	+	+	-	+	+
Methanol	1%	-	+	+	-	+	+
Heparin	1%	-	+	+	-	+	+
Citrate	3.2%	-	+	+	-	+	+
Biological analytes							
Albumin	2 mg/ml	-	+	+	-	+	+
Glucose	2 mg/ml	-	+	+	-	+	+
Bilirubin	2 mg/ml	-	+	+	-	+	+
Hemoglobin	2 mg/ml	-	+	+	-	+	+

REFERENCES

1. Wolters G, Kuijpers L, Kacaki J, Schuurs AHWM. Solid-phase enzyme-immunoassay for detection of hepatitis B surface antigen. *J Clin Pathol*(1976), 29 (10): 873-879.
2. David GS, Present W, Martinis J, Wang R, Bartholomew R, Desmond W, Sevier ED. Monoclonal antibodies in the detection of hepatitis infection. *Med Lab Sci* (1981)38: 341-348.
3. De Meyer S, Gong ZJ, Suwandhi W, van Pelt J, Soumillion A, Yap SH. Organ and species specificity of hepatitis B virus (HBV) infection: a review of literature with a special reference to preferential attachment of HBV to human hepatocytes. *J Viral Hepat* (1997), 4:145-153.
4. Cho EW, Park JH, Yoo OJ, Kim KL. Translocation and accumulation of exogenous hepatitis B virus preS surface proteins in the cell nucleus. *J Cell Sci*(2001). 114 (Pt 6):1115-1123.
5. Cooper A, Paran N, Shaul Y. The earliest steps in hepatitis B virus infection. *Biochim Biophys Acta* (2003). 1614: 89-96.

INDEX OF SYMBOLS

	Do not reuse		Batch code
	In vitro diagnostic medical device		Use by
	Temperature limitation		Contains sufficient for < n > tests
	Caution		Catalog number
	Manufacturer		Consult instructions for use
	Authorized representative in the European community		CE Mark

MANUFACTURER CONTACT INFORMATION

Atlas Link (Beijing) Technology Co., Ltd
Room 811 Zeyang Plaza. No.166 Fushi Road
Shijingshan Dist. Beijing
100043, China
website: <https://www.invitro-test.com>

Tel: +86-10-88909113



CIRIANO GLOBAL S.L (UNILATEX GROUP)
CIF: B50927532
C/Blancas 4-6
50001 Zaragoza, SPAIN

Tel: +34-976228974

Doc No. ATARHBV211

REVISION: May, 2014