

One Step Cancer CEA Serum Rapid Test Strip

INTRODUCTION

One Step Strip Style CEA Test is a rapid direct binding test for the detection of Carcinoma Embryonic Antigen (CEA) in serum as an aid in the diagnosis of colorectal, breast, lung and pancreas cancers. The test is based on the principle of sandwich immunoassay for determination of CEA in serum. Monoclonal and polyclonal antibodies are employed to identify CEA specifically. This one step test is very sensitive and only takes about 10-20 minutes. The sensitivity of the test can reach to 4 µg/L.


SPECIMEN COLLECTION

For serum, collect blood into a container without anticoagulant. Allow the blood to clot and separate the serum from the clot. Use the serum for testing.

If the specimen cannot be tested on the day of collection, store the serum specimen in a refrigerator or freezer. Bring the specimens to room temperature before testing. Do not freeze and thaw the specimen repeatedly.


TEST PROCEDURE

1. Following the illustration, dip the strip with arrow pointing down into the vessel of serum for about 10 seconds. Do not immerse past the MAX line.
2. Wait for 10 minutes and read results. Do not read results after 30 minutes.



INTERPRETATION OF RESULTS

- **Negative:** Only one colored band appears on the control region. No apparent band on the test region.
- **Positive:** Distinct color bands appear on the control and test regions. Both test line and control line indicate a positive CEA result. Color intensity of the test bands may vary.
- **Invalid:** A total absence of color in both regions is an indication of procedure error and/or that the test reagent has deteriorated.



STORAGE AND STABILITY

The test kits can be stored at room temperature (2 to 30°C) in the sealed pouch to the date of expiration. The test kits should be kept away from direct sunlight, moisture and heat.

PRECAUTION

1. For in vitro diagnostic use only.
2. Do not use test kit beyond the expiry date.
3. The test device should not be reused.

MANUFACTURER

ATLAS LINK TECHNOLOGY CO., LTD

Road Xing Min, Guan South Industry Zone, 065500 Langfang City, Hebei Province, CHINA

Web: <https://www.invitro-test.com>